



Diagnostics – Anesthesia
Diagnostik – Anästhesie
Diagnóstico – Anestesia
Diagnostic – Anesthésie
Diagnostica – Anestesia



Babinsky
 01-100-22
 22 cm - 8 ¾"



Rabiner
 01-101-23
 23 cm - 9 ¼"



Queen's Hammer
 01-108-35
 35 cm - 13 ¾"
 Plastic handle



Wartenberg
01-102-19
19 cm - 7 ½"



Buck
01-103-18
18 cm - 7"



Traube
01-104-18
18 cm - 6 ¼"



Taylor
01-105-18
18 cm - 7"
01-105-20
20 cm - 8"



Berliner
01-106-20
20 cm - 8"



Déjerine
01-107-20
20 cm - 8"

Set 1
01-108-00

Ophthalmoscope set
Ophthalmoskop-Besteck
Ophtalmoscopio
Trousse à ophtalmoscopie
Oftalmoscopio



Set 2
01-109-00

Universal set
Universal-Besteck
Estuche universal
Trousse universelle
Astuccio universale



Set 3
01-110-00

Combination set
Kombinations-Besteck
Estuche combinado
Trousse combinée
Astuccio combinato



Set 4
01-111-00

Otoscope set
Otoskop-Besteck
Otoscopio
Trousse à otoscopie
Otoscopio





Pinard
01-112-18
18 cm - 7"
Wood
Holz
Madera
Bois
Legno



Pinard
01-113-15
15 cm - 6"
Aluminium
Aluminium
Aluminio
Aluminium
Alluminio



Duplex
01-114-04 4.7cm - 1 3/4"
for adults
für Erwachsene
para adultos
pour adultes
per adulti

Baby-Duplex
01-115-15 4.7cm - 1 3/4"
for children
für Kinder
para niños
pour enfants
per bambini



Anestophon
01-116-04 4.7cm - 1 3/4"
for nurses, simple
für Schwestern, einfach
para enfermeras, simple
pour infirmières, simple
per infermiere, semplice





Ziegler

01-117-09 9 cm - 3 ½" Ø

Head mirror with black plastic head band
Stirnspiegel mit schwarzem Plastikband
Espejo frontal con banda de plastica negra
Miroir frontal avec ruban en plastique noir
Specchio frontale con fascia di plastica nera

01-119-09 9 cm - 3 ½" Ø

Head mirror with white plastic head band
Stirnspiegel mit weissem Plastikband
Espejo frontal con banda de plastica blanca
Miroir frontal avec ruban en plastique blanc
Specchio frontale con fascia di plastica bianca

Ziegler

01-118-09 9 cm - 3 ½" Ø

Head mirror only, with rubber coating
Reflexspiegel allein, mit Gummifassung
Espejo frontal solo, con guarnición de goma
Miroir frontal seul, avec monture en caoutchouc
Solo specchio frontale, con guarnizione di gomma

01-120-00

Black plastic head band
Plastikband, schwarz
Banda de plastica, negra
Ruban en plastique noir
Fascia di plastica, nera

01-121-00

White plastic head band
Plastikband, weiss
Banda de plastica, blanca
Ruban en plastique blanc
Fascia di plastica, bianca



Cold-light head lamp
Kaltlicht-Stirnlampe
Lámpara frontal de luz fría
Eclaireur frontal de lumière froide
Lampada frontale a luce fredda

Binner

01-122-00

Head light, cold-light, adjustable focus,
working distance 300 mm, field of illumination 70 – 120 mm Ø

Kaltlicht-Stirnlampe mit verstellbarer Fokuserlinse,
Arbeitsabstand 300mm, Leuchtfeldgröße 70 – 120 mm Ø

Lámpara frontal, luz fría, focus ajustable,
distancia de trabajo 300mm, campo de iluminación 70 – 120mm Ø

Eclaireur frontal, lumière froide, focus ajustable,
distance de travail 300 mm, champ d'illumination 70 – 120mm Ø

Lampada frontale, luce fredda, con lente a focale variabile,
distanza di lavoro 300 mm, campo d'illuminazione 70 – 120mm Ø



Clar

01-123-00 10 cm Ø, 6 V

Head light complete
Stirnlampe komplett
Lámpara frontal, completa
Eclaireur frontal complet
Lampada frontale, completa

01-124-00 10 cm Ø, 6 V

Head light with cable, frame and mirror, without head band
Stirnlampe mit Kabel, Fassung und Spiegel, ohne Stirnreifen
Lámpara frontal, con cable, armadura y espejo, sin banda
Eclaireur frontal avec cable, douille et miroir, sans ruban
Lampada frontale con cavo, supporto e specchio, senza fascia



Mini-Clar

01-125-00 5 cm Ø, 6 V

complete
komplett
completa
complet
completa

Clar

01-126-05 5 cm - 2" Ø

01-126-10 10 cm - 2" Ø

Spiegel aus Glas, ohne Fassung
Head mirror glass, without frame
Miroir frontal en vitre, sans douille
Espejo frontal de cristal, sin armadura
Vetro specchio, senza supporto

Clar

01-127-00

Spezial-Gelenkverbindung
Special socket joint
Assemblage articulé spécial
Conexión
Giunto di connessione speciale

Clar

01-128-00

Stirnreifen allein
Head band only
Ruban seul
Banda sola
Solo fascia frontale

01-129-00 6 V

Glühbirne
Bulb
Ampoule
Bombilla
Lampadina



Transformer – stepless adjustable
Transformator – stufenlos einstellbar
Transformador – regulación continua
Transformateur – à réglage continu
Trasformatore – regolabile in continuità

01-130-00

Input 220 V/ 50 Hz ~
Output 0–10 V/ 3 W =

Binocular

01-131-00

2.5 times enlargement, 340 mm working distance, field of vision 90 mm Ø, complete.
Pupillary distance adjustable from 48 mm–77 mm. Individual adjustable angle of inclination.
Weight of the optic 42g.

2.5 fache Vergrößerung, 340 mm Arbeitsabstand, Sichtfeld > 90 mm Ø, komplett.
Pupillenabstand einstellbar von 48 mm–77 mm. Neigungswinkel individuell justierbar.
Gewicht der Optik 42g.

2.5 veces aumento, distancia de trabajo 340 mm, campo de visión >90 mm Ø, completo.
Distancia interpupilar ajustable de 48 mm–77 mm. Angulo de inclinación ajustable individual.
Peso de la optica 42g.

2.5 fois grossissement, distance de travail 340 mm, champ de vision >90 mm Ø, complèt.
Ecart interpupillaire ajustable de 48 mm–77 mm. Angle d'inclinaison ajustable individuellement.
Poids de l'optique 42g.

Ingrandimento: 2.5 volte, distanza di lavoro 340 mm, campo di visione >90 mm Ø, completo.
Distanza interpupillare regolabile da 48 mm–77 mm. Angolo di inclinazione regolabile individualmente.
Peso dell'ottica 42g.



01-132-00

Kopfband
Head band
Ruban
Banda
Fascia

Sphygmomanometers
Blutdruckmessgeräte
Esfigmomanetros
Tensiomètres
Sfigmomanometri

Manometer Fig. 1

01-133-00

Calibrated
Geeicht
Calibado
Calibré
Calibrato



Manometer Fig. 2

01-134-00

non calibrated
ungeeicht
no calibrado
non calibré
non calibrato



Manometer Fig. 3

01-135-00

Mercury blood pressure manometer, non calibrated
Quecksilber-Blutdruckmessapparat, ungeeicht
Manómetro de mercurio, no calibrado
Tensiomètre à mercure, non calibré
Manometro a mercurio, non calibrato





Collin
 Fig. 1

01-136-12

12.5 cm - 5"

for children, kg-scale
 für Kinder, kg-Einteilung
 para niños, graduación en kg
 pour enfants, échelle en kg
 per bambini, graduazione in kg

Fig. 2

01-137-12

12.5 cm - 5"

for adults, kg-scale
 für Erwachsene, kg-Einteilung
 para adultos, graduación en kg
 pour adultes, échelle en kg
 per adulti, graduazione in kg

Fig. 3

01-138-12

12.5 cm - 5"

for adults, lbs-scale
 für Erwachsene, lbs-Einteilung
 para adultos, graduación en lbs
 pour adultes, échelle en lbs
 per adulti, graduazione in lbs



Bertillon

01-139-32

32 cm - 12 1/2"



Moeltgen

01-140-24

24 cm 9 1/2"



Magda
Fig. 1

01-141-10	10 cm - 4"
01-141-15	15 cm - 6"
01-141-20	20 cm - 8"
01-141-30	30 cm - 12"
01-141-50	50 cm - 20"



Magda
Fig. 2

01-142-50	50 cm - 20"
X-ray ruler Röntgen-Messstab Regla para rayos X Règle rayons X Misuratore per raggi X	



Vernier

01-143-21	01-144-21	21 cm - 8 1/4"
Aluminium	Steel	
Aluminium	Stahl	
Aluminio	Acero	
Aluminium	Acier	
Alluminio	Acciaio	



Castroviejo

01-145-08
8 cm - 3 1/4"



Gauge

01-146-15	01-147-15	150 cm - 60"
Steel	Wax cloth	
Stahl	Wachs	
Acero	Cerada	
Acier	Toile cirée	
Acciaio	Cera	



Welzer

01-152-16 16 cm - 6 ¼"

Wood
Holz
Madera
Bois
Legno



Wavy

01-153-14 14.5 cm - 5 ¾"



Kroenig

01-154-15 15 cm - 6"



Mayo

01-155-17 17 cm - 6 ¾"



Bruenings

01-156-19
19 cm - 7 ½"



Frenzel

01-157-16
16.5 cm - 6 ½"



Lack

01-130-13	13 mm
01-130-19	19 mm
01-130-22	22 mm



Wieder

01-158-15
15 cm - 6"



01-158-14
14 cm - 5 ½"





Bosworth
01-159-14
14.5 cm - 5 ¾"



Tobold
01-160-14
14.5 cm - 5 ¾"



Hartmann
01-161-15
15 cm - 6"



Andrew
01-162-11
11 cm - 4 ¾"



Young
01-163-17
17 cm - 6 ¾"



01-163-00
1 Pair rubber-inserts
1 Paar Ersatz Gummi-Einsätze
1 Juego de insertos de goma
1 Paire des inserts en caoutchouc
1 Paio d'inserti di gomma



Collin
01-164-16 16 cm - 6 ¼"
01-164-17 17 cm - 6 ¾"
01-164-19 19 cm - 7 ½"

Cheek retractors
Wangenhalter
Separador de mejillas
Ecarteur des joues
Divaricatori per guancia



Sternberg

01-165-12	12 cm - 4 ¾"
01-165-14	14 cm - 5 ½"
01-165-16	16 cm - 6 ¼"

Catheter forceps
Katheterzangen
Pinzas para catéteres
Pincers pour cathéters
Pinze per catetere



Magill

01-166-17

17 cm - 6"
Infant



01-166-20

20 cm - 8"
Child



01-166-25

25 cm - 9"
Adult



The English Macintosh is similar to the conventional Macintosh except that the flange is curved and lower at the handle end. Two studies have proved that this blade provided better results than other commonly used blades.



Fig. 00



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5

Code	Description	Fig.	Size
01-281-00	Premature	00	68 mm
01-281-0	Neonate	0	80 mm
01-281-1	Infant	1	92 mm
01-281-2	Child	2	105 mm
01-281-3	Medium Adult	3	135 mm
01-281-4	Large Adult	4	155 mm
01-281-5	X-Large	5	175 mm

The Miller is one of the most popular blades. The tongue is straight with a slight upward curve near the tip. In cross-section the flange, web, and tongue form a C with the top fattened. Placing the bulb on the right side allows it to be protected by the flange. Several modifications have been described. In comparing the Macintosh and Miller blades, it was found that force, head extension and cervical spine movement were less with the Miller.



Fig. 00



Fig. 2



Fig. 0



Fig. 3



Fig. 1



Fig. 4

Code	Description	Fig.	Size
01-282-00	Premature	00	65 mm
01-282-0	Neonate	0	78 mm
01-282-1	Infant	1	102 mm
01-282-2	Child	2	154 mm
01-282-3	Medium Adult	3	195 mm
01-282-4	Large Adult	4	205 mm

Specially designed and constructed with high quality stainless steel AISI 304.
High quality German optical fibers loaded with 4.5mm bundle for best optimum light transmission.

Exceptionally high light intensity with LED illumination (8000-10000 LUX).

With this integrated FO light conductor, there are no extra cavities or light path to become contaminated.

Smooth and low profile blade design which gives no chance of damaging teeth and upper part of mouth.

Integrated light guide blade design for quick cleaning and sterilization without extra effort for connecting screw.

Compatible with all handles of ISO 7376 (Green System).

Autoclaving at 134 °C / 5 min. maximum light output with 3000-4000 cycles.



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5



Code	Description	Fig.	Size
01-283-0	Neonate	0	80 mm
01-283-1	Infant	1	93 mm
01-283-2	Child	2	114 mm
01-283-3	Medium Adult	3	135 mm
01-283-4	Large Adult	4	155 mm
01-283-05	X-Large	5	175 mm

Specially designed and constructed with high quality stainless steel AISI 304.
 High quality German optical fibers loaded with 4.5mm bundle for best optimum light transmission.

Exceptionally high light intensity with LED illumination (8000-10000 LUX).

With this integrated FO light conductor, there are no extra cavities or light path to become contaminated.

Smooth and low profile blade design which gives no chance of damaging teeth and upper part of mouth.

Integrated light guide blade design for quick cleaning and sterilization without extra effort for connecting screw.

Compatible with all handles of ISO 7376 (Green System).

Autoclaving at 134 °C / 5 min. maximum light output with 3000-4000 cycles.



Fig. 00



Fig. 2



Fig. 0



Fig. 3



Fig. 1



Fig. 4

Code	Description	Fig.	Size
01-284-00	Premature	00	65 mm
01-284-0	Neonate	0	78 mm
01-284-1	Infant	1	102 mm
01-284-2	Child	2	154 mm
01-284-3	Medium Adult	3	195 mm
01-284-4	Large Adult	4	205 mm

These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered.

It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis.




Fig. 1



Fig. 2



Fig. 3

Code	Description	Fig.	Size
01-285-2	Child	2	105 mm
01-285-3	Medium Adult	3	130 mm
01-285-4	Large Adult	4	155 mm

These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered.

It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. However, it will not always improve the view and may worsen it. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis. In some patients in whom the view is satisfactory with the Macintosh blade, use of the Flexible tip blade may make the view worse.



Fig. 2



Fig. 3



Fig. 4

Code	Description	Fig.	Size
01-286-2	Child	2	108 mm
01-286-3	Medium Adult	3	130 mm
01-286-4	Large Adult	4	155 mm

Diamond knurled barrel gives strong grip to the user.
LED brighter illumination with 60,000 LUX.
LED operating time over 50,000 hrs.
2 hrs operating time at maximum illumination.
Viewing angle at 130 with color temperature color temperature 6500K.
Conforms to all Green System blades as per ISO 7376 standard.
Simple change of batteries at the base of the handle.
Crafted from high quality chrome plated brass.
Battery usage is minimized due to the low power consumption.
Fully autoclavable @ 134 C for 5 minutes.

Battery Handles Fiber Optic Xenon, Stainless Steel

Code	Description	Battery	Size
01-287-01	Infant	1xLithium 123 Size	122 mm
01-287-02	Stubby	2xAA Size	160 mm
01-287-03	Pediatric	2xAA Size	125 mm
01-287-04	Medium	2xC Size	160 mm

Battery Handles Fiber Optic Xenon, Brass Chrome Plated

Code	Description	Battery	Size
01-288-01	Infant	1xLithium 123 Size	122 mm
01-288-02	Stubby	2xAA Size	160 mm
01-288-03	Pediatric	2xAA Size	125 mm
01-288-04	Medium	2xC Size	160 mm



Battery Handles Fiber Optic LED, Stainless Steel

Code	Description	Battery	Size
01-289-01	Infant	1xLithium 123 Size	122 mm
01-289-02	Stubby	2xAA Size	160 mm
01-289-03	Pediatric	2xAA Size	125 mm
01-289-04	Medium	2xC Size	160 mm

Battery Handles Fiber Optic LED, Brass Chrome Plated

Code	Description	Battery	Size
01-290-01	Infant	1xLithium 123 Size	122 mm
01-290-02	Stubby	2xAA Size	160 mm
01-290-03	Pediatric	2xAA Size	125 mm
01-290-04	Medium	2xC Size	160 mm



Macintosh Flaxiblade Pro Blades, Fiber Optic, German Profile
 Macintosh Flaxiblade Pro Klingen, Faseroptik, Deutsches Profil
 Macintosh Flaxiblade Pro Blades, Fibra Óptica, Perfil Alemán
 Macintosh Flaxiblade Pro Lames, fibre optique, profil allemand
 Macintosh Flaxiblade Pro Lame, Fibra ottica, Profilo tedesco



Battery Handles Fiber Optic LED, Stainless Steel

Code	Description	Battery	Size
01-291-01	Infant	1xLithium 123 Size	122 mm
01-291-02	Stubby	2xAA Size	160 mm
01-291-03	Pediatric	2xAA Size	125 mm
01-291-04	Medium	2xC Size	160 mm



LED Lamp



Halogen Lamp

Replacement Lamps

Code	Type	Voltage
01-292-01	Xenon	2.5 V
01-292-02	LED	3.0 V
01-292-03	Halogen	2.5 V
01-292-04	Halogen	3.5 V

Replacement Fiber Optic Light Guides, English Profile (4 mm)

Code	Description	For Blades
01-296-00	Miller Blade	00
01-296-0	Miller Blade	0
01-296-1	Miller Blade	1
01-296-2	Miller Blade	2
01-296-3	Miller Blade	3
01-296-4	Miller Blade	4
01-297-00	Macintosh Blade	00
01-297-0	Macintosh Blade	0
01-297-1	Macintosh Blade	1
01-297-2	Macintosh Blade	2
01-297-3	Macintosh Blade	3
01-297-4	Macintosh Blade	4



Macintosh Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
01-298-01	Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-298-02	Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-298-03	Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Miller Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
01-299-01	Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-299-02	Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-299-03	Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Flexible Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
01-300-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

Flexible pro Fiber Optic Xenon Laryngoscope Sets, German Profile

Code	Description
01-301-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Macintosh Fiber Optic LED Laryngoscope Sets, English Profile

Code	Description
01-302-01	Mac F.O LED Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-302-02	Mac F.O LED Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-302-03	Mac F.O LED Laryngoscope Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Miller Fiber Optic LED Laryngoscope Sets, English Profile

Code	Description
01-303-01	Miller F.O LED Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-303-02	Miller F.O LED Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-303-03	Miller F.O LED Laryngoscope Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Flexible Fiber Optic Led Laryngoscope Sets, English Profile

Code	Description
01-304-01	MacFlex F.O LED Laryngoscope Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

Flexible pro LED Laryngoscope Sets, German Profile

Code	Description
01-305-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

The English Macintosh is similar to the conventional Macintosh except that the flange is curved and lower at the handle end. Two studies have proved that this blade provided better results than other commonly used blades.



Code	Description	Fig.	Size
01-306-00	Premature	00	68 mm
01-306-0	Neonate	0	80 mm
01-306-1	Infant	1	95 mm
01-306-2	Child	2	105 mm
01-306-3	Medium Adult	3	135 mm
01-306-35	Ext. Medium Adult	3.5	145 mm
01-306-4	Large Adult	4	155 mm
01-306-5	X-Large	5	175 mm

The Miller is one of the most popular blades. The tongue is straight with a slight upward curve near the tip. In cross-section the flange, web, and tongue form a C with the top fattened. Placing the bulb on the right side allows it to be protected by the flange. Several modifications have been described. In comparing the Macintosh and Miller blades, it was found that force head extension and cervical spine movement were less with the Miller.



Fig. 00



Fig. 2



Fig. 0



Fig. 3



Fig. 1



Fig. 4

Code	Description	Fig.	Size
01-307-00	Premature	00	65 mm
01-307-0	Neonate	0	78 mm
01-307-1	Infant	1	102 mm
01-307-2	Child	2	154 mm
01-307-3	Medium Adult	3	195 mm
01-307-4	Large Adult	4	205 mm

These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered. It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis.



Fig. 2



Fig. 3



Fig. 4

Code	Description	Fig.1	Size
01-308-2	Child	2	105 mm
01-308-3	Medium Adult	3	130 mm
01-308-4	Large Adult	4	155 mm

Battery Handles & Replacement Lamps for Conventional Laryngoscope
 Batteriegriffe & Ersatzlampen für herkömmliches Laryngoskop
 Manijas de baterías y lámparas de repuesto para laringoscopia convencional
 Poignées de batterie et lampes de rechange pour le laryngoscope conventionnel
 Batteria Maniglie e lampade di ricambio per laringoscopia convenzionale

Diamond knurled barrel gives strong grip to the user.
 Conforms to all conventional blades as per ISO 7376 standard.
 Simple change of batteries at the base of the handle.
 Crafted from high quality chrome plated brass.
 Economical solution for intubation as compare to green system.
 Fully autoclavable @ 134 °C for 5 minutes.




Battery Handles, Stainless Steel

Code	Description	Battery	Size
01-309-01	Infant	1xLithium 123 Size	110 mm
01-309-02	Stubby	2xAA Size	115 mm
01-309-03	Pediatric	2xAA Size	145 mm
01-309-04	Medium	2xC Size	150 mm

Battery Handles, Brass Chrome Plated

Code	Description	Battery	Size
01-310-01	Infant	1xLithium 123 Size	110 mm
01-310-02	Stubby	2xAA Size	115 mm
01-310-03	Pediatric	2xAA Size	145 mm
01-310-04	Medium	2xC Size	150 mm

Replacement Lamps



Code	Description	Type	Voltage
01-311-00	Small	Xenon	2.5 V
01-312-00	Adult	LED	3.0 V
01-313-00	Small	Halogen	2.5 V
01-314-00	Adult	Halogen	3.5 V



Macintosh Conventional Laryngoscope Sets, English Profile

Code	Description
01-315-01	Mac Conventional Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-315-02	Mac Conventional Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-315-03	Mac Conventional Laryngoscope Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Miller Conventional Laryngoscope Sets, English Profile

Code	Description
01-316-01	Miller Conventional Laryngoscope Sets Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp in Plastic Case.
01-316-02	Miller Conventional Laryngoscope Sets Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
01-316-03	Miller Conventional Laryngoscope Sets Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Jackson

Jackson with light carrier, stainless steel
Jackson mit Lichtträger, rostfrei
Jackson con porta fibras, inoxidable
Jackson avec porte de lumière, inoxydable
Jackson con porta luce, inossidabile

01-188-12

x = 12 cm
for babys
für Säuglinge
para lactantes
pour nourissons
per lattanti

01-188-17

x = 17 cm
for children
für Kinder
para niños
pour enfants
per bambini

01-188-18

x = 18 cm
for adults
für Erwachsene
para adultos
pour adultes
per adulti

Light carrier, only, for Jackson
Lichtträger, allein, für Jackson
Porta fibras, sola, para Jackson
Porte de lumière, seul, pour Jackson
Porta luce, sola, per Jackson

01-189-12

x = 12 cm
for babys
für Säuglinge
para lactantes
pour nourissons
per lattanti

01-189-17

x = 17 cm
for children
für Kinder
para niños
pour enfants
per bambini

01-189-18

x = 18 cm
for adults
für Erwachsene
para adultos
pour adultes
per adulti



Holinger

01-190-17

x = 17.1 cm

01-190-18

x = 18 cm

Holinger with light carrier, stainless steel
Holinger mit Lichtträger, rostfrei
Holinger con porta fibras, inoxidable
Holinger avec porte de lumière, inoxydable
Holinger con porta luce, inossidabile

01-191-17

x = 17.1 cm

01-191-18

x = 18 cm

Light carrier, only, for Holinger
Lichtträger, allein, für Holinger
Porta fibras, sola, para Holinger
Porte de lumière, seul, pour Holinger
Porta luce, solo, per Holinger



Inka
Fig.1

01-148-01 x = 13.5 cm

Stainless steel
Rostfreier Edelstahl
Acero inoxidable
Acier inoxydable
Acciaio inossidabile



Inka
Fig.2

01-149-02 x = 15.2 cm

Stainless steel
Rostfreier Edelstahl
Acero inoxidable
Acier inoxydable
Acciaio inossidabile



Inka
Fig.3

01-150-03 x = 18.2 cm

Stainless steel
Rostfreier Edelstahl
Acero inoxidable
Acier inoxydable
Acciaio inossidabile



Inka
Fig.4

01-151-04 x = 17.2 cm

Stainless steel
Rostfreier Edelstahl
Acero inoxidable
Acier inoxydable
Acciaio inossidabile



Inka

01-192-05 x = 17.2 cm

Fig. 5

Stainless steel
 Rostfreier Edelstahl
 Acero inoxidable
 Acier inoxydable
 Acciaio inossidabile



Inka

01-193-06 x = 18.2 cm

Fig. 6

Stainless steel
 Rostfreier Edelstahl
 Acero inoxidable
 Acier inoxydable
 Acciaio inossidabile



Inka

01-194-07 x = 22.2 cm

Fig. 7

Rostfreier Edelstahl
 Stainless steel
 Acier inoxydable
 Acero inoxidable
 Acciaio inossidabile



Jason


01-195-34 34 cm - 13"

Chest-support with support-rod for adults
 Bruststütze mit Stützstange für Erwachsene
 Soporte de pecho con vara para adultos
 Appui-poitrine avec tige pour adultes
 Supporto per petto con asta per adulti



01-196-24 24 cm - 9"

Support-rod for children
 Stützstange für Kinder
 Vara para niños
 Tige pour enfants
 Asta per bambini



Lupus
Fig. 1

01-197-00

for laryngoscopes
für Laryngoskope
para laringoscopios
pour laryngoscopes
per laringoscopi




Fig. 2

01-198-00

Suction tube for laryngoscopes
Rauchabsaugrohr für Laryngoskope
Tubo de aspiración para laringoscopios
Tube d'aspiration pour laryngoscopes
Tubo d'aspirazione per laringoscopi



Fig. 3

01-199-00

for laryngoscopes
für Laryngoskope
para laringoscopios
pour laryngoscopes
per laringoscopi



Fig. 4

01-200-00

Injection drain tube for laryngoscopes
Injektionskanüle für Laryngoskope
Canula de inyección para laringoscopios
Canule à injection pour laryngoscopes
Cannula d'iniezione per laringoscopi




Fig. 5

01-201-00

for laryngoscope
für Laryngoskop
para laringoscopia
pour laryngoscope
per laringoscopia