





Sterilization
Sterilisation
Esterilización
Stérilisation
Sterilizzazione

comed
operate better with us

21

Packing		Requirements according to:						
Type	Description	Symbolic Delivery and Packaging System	DIN EN 285	DIN 58946-1				
				Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure
	Textiles in sterilization container with filter inside lid and bottom.		+	+	+	+	+	+
			+	0	0	+	+	+
	Textiles in sterilization container with filter inside lid, and not perforated bottom.		+	-	-	+	+	+
			+	-	-	+	+	0
	Textiles in sterilization container with filter inside bottom and not perforated lid.		+	0	+	+	+	+
			+	-	0	+	+	+
	Instruments in wire baskets inside sterilisation container with filter inside lid and bottom.		+	+2)	+	+	+	+
	Instruments in wire baskets in sterilization container with filter inside lid and not perforated bottom.		+	-	+2)	+2)	+2)	+3)
	Instruments in wire baskets inside sterilization container with filter inside bottom.		+	0	+	+	+	+

Packing		Requirements according to:						
Type	Description	Symbolic Delivery and Packaging System	DIN EN 285	Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure
	Textiles in sterilization container with valve inside lid and bottom.		+	-	0	+	+	0
			+	-	-	+	+	0
	Textiles in sterilization container with valve inside lid and no valve in the bottom.		+	-	-	+	+	-
			+	-	-	+	+	-
	Instruments in wire baskets inside sterilization container with valve inside lid and bottom.		+	-	+2)	+2)	+2)	+2)
	Instruments in wire baskets inside sterilization container with valve inside lid, no valve in the bottom.		+	-	+2)	+2)	+2)	+2)
Legend: + applicable, – not applicable, 0 applicable under certain conditions								
1). To avoid service, operation and/or loading errors, see Table 1 for some recommendations regarding the correct or incorrect combination of sterilization containers with standard steam sterilisation procedures.								
2). Only sterilisation containers from which the condensate evaporation or flow off is guaranteed are suitable.								
3). The gravity procedure without a drying device does not provide sufficient re-drying.								



Standard Model



Safety Model

GB

Choose between our standard and our safety model. The safety lid protects against contamination while transport and stocking. Filter- and safety lid can only separate by opening container.

DE

Wählen Sie zwischen unserem Standard- und unserem Sicherheitsmodell. Der Sicherheitsdeckel schützt vor Kontamination bei Transport und Lagerung. Filter- und Sicherheitsdeckel lassen sich nur durch Öffnen des Behälters trennen.

ES

Elija entre nuestro estándar y nuestro modelo de seguridad. La tapa de seguridad protege contra la contaminación durante el transporte y almacenamiento. El filtro y la tapa de seguridad solo pueden separarse abriendo el recipiente.

FR

Choisissez entre notre modèle standard et notre modèle de sécurité. Le couvercle de sécurité protège contre la contamination pendant le transport et le stockage. Le filtre et le couvercle de sécurité ne peuvent se séparer qu'en ouvrant le récipient.

IT

Scegli tra il nostro modello standard e il nostro modello di sicurezza. Il coperchio di sicurezza protegge dalla contaminazione durante il trasporto e lo stoccaggio. Il filtro e il coperchio di sicurezza possono essere separati solo aprendo il contenitore.



Lid Perforated, Bottom Non Perforated



Lid Perforated, Bottom Perforated



GB

The locking flaps are movable 180° on a pivot and can be secured in any position.

DE

Die Verschlussklappen sind um 180° auf einem Drehpunkt beweglich und können in jeder Position arretiert werden.

ES

Las aletas de bloqueo se pueden mover 180° sobre un pivote y se pueden asegurar en cualquier posición.

FR

Les volets de verrouillage sont mobiles à 180° sur un pivot et peuvent être bloqués dans n'importe quelle position.

IT

Le alette di chiusura sono mobili di 180° su un perno e possono essere fissate in qualsiasi posizione.



GB

The filter retainer is locked by pressing down in center of the filter. Opening the filter retainer by pressing.

DE

Der Filterhalter wird durch Herunterdrücken in der Mitte des Filters verriegelt. Öffnen des Filterhalters durch Drücken.

ES

El retenedor del filtro se bloquea presionando hacia abajo en el centro del filtro. Apertura del portafiltros presionando.

FR

Le dispositif de retenue du filtre est verrouillé en appuyant au centre du filtre. Ouverture du porte-filtre par pression.

IT

Il fermo del filtro si blocca premendo al centro del filtro. Apertura del portafiltro premendo.

1/1 Size Container - Standard Model



Non perforated
Bottom



Perforated
Bottom

	21-444-01
	21-444-02
	21-444-03
	21-444-04
	21-444-05
	21-444-06

Weight: 2634 gm

	21-445-01
	21-445-02
	21-445-03
	21-445-04
	21-445-05
	21-445-06

Weight: 3156 gm



Outside: 580 x 280 x 100 mm
 Inside: 550 x 265 x 84 mm

	21-446-01
	21-446-02
	21-446-03
	21-446-04
	21-446-05
	21-446-06

Weight: 3125.5 gm

	21-447-01
	21-447-02
	21-447-03
	21-447-04
	21-447-05
	21-447-06

Weight: 3244 gm



Outside: 580 x 280 x 135 mm
 Inside: 550 x 265 x 115 mm

	21-444-01
	21-444-02
	21-444-03
	21-444-04
	21-444-05
	21-444-06

Weight: 2634 gm

	21-445-01
	21-445-02
	21-445-03
	21-445-04
	21-445-05
	21-445-06

Weight: 3156 gm



Outside: 580 x 280 x 150 mm
 Inside: 550 x 265 x 130 mm

	21-448-01
	21-448-02
	21-448-03
	21-448-04
	21-448-05
	21-448-06

Weight: 3219 gm

	21-449-01
	21-449-02
	21-449-03
	21-449-04
	21-449-05
	21-449-06

Weight: 3382 gm



Outside: 580 x 280 x 200 mm
 Inside: 550 x 265 x 180 mm

	21-450-01
	21-450-02
	21-450-03
	21-450-04
	21-450-05
	21-450-06

Weight: 3429.5 gm

	21-451-01
	21-451-02
	21-451-03
	21-451-04
	21-451-05
	21-451-06

Weight: 3738.5 gm



Outside: 580 x 280 x 260 mm
 Inside: 550 x 265 x 240 mm

1/1 Size Container - Safety Model



Non perforated
Bottom



Perforated
Bottom

	21-454-01
	21-454-02
	21-454-03
	21-454-04
	21-454-05
	21-454-06

Weight: 3486 gm

	21-455-01
	21-455-02
	21-455-03
	21-455-04
	21-455-05
	21-455-06

Weight: 4007.5 gm



Outside: 580 x 280 x 100 mm
Inside: 550 x 265 x 84 mm

	21-456-01
	21-456-02
	21-456-03
	21-456-04
	21-456-05
	21-456-06

Weight: 3981.5 gm

	21-457-01
	21-457-02
	21-457-03
	21-457-04
	21-457-05
	21-457-06

Weight: 4095.5 gm



Outside: 580 x 280 x 135 mm
Inside: 550 x 265 x 115 mm

	21-458-01
	21-458-02
	21-458-03
	21-458-04
	21-458-05
	21-458-06

Weight: 4070.5 gm

	21-459-01
	21-459-02
	21-459-03
	21-459-04
	21-459-05
	21-459-06

Weight: 4233 gm



Outside: 580 x 280 x 150 mm
Inside: 550 x 265 x 130 mm

	21-460-01
	21-460-02
	21-460-03
	21-460-04
	21-460-05
	21-460-06

Weight: 4280 gm

	21-461-01
	21-461-02
	21-461-03
	21-461-04
	21-461-05
	21-461-06

Weight: 4589 gm



Outside: 580 x 280 x 200 mm
Inside: 550 x 265 x 180 mm

	21-462-01
	21-462-02
	21-462-03
	21-462-04
	21-462-05
	21-462-06

Weight: 4770.5 gm

	21-463-01
	21-463-02
	21-463-03
	21-463-04
	21-463-05
	21-463-06

Weight: 5140 gm



Outside: 580 x 280 x 260 mm
Inside: 550 x 265 x 240 mm



Non perforated
Bottom

	Outside	Inside	Weight
21-464-01	580 x 280 x 100 mm	550 x 262 x 84 mm	1923 gm
21-464-02	580 x 280 x 135 mm	550 x 262 x 117 mm	2019 gm
21-464-03	580 x 280 x 150 mm	550 x 262 x 132 mm	2108 gm
21-464-04	580 x 280 x 200 mm	550 x 262 x 182 mm	2317 gm
21-464-05	580 x 280 x 260 mm	550 x 262 x 244 mm	2808 gm

1/1 Size Container - Bottoms and Lids



Perforated
Bottom

	Outside	Inside	Weight
21-465-01	580 x 280 x 100 mm	550 x 262 x 84 mm	2044 gm
21-465-02	580 x 280 x 135 mm	550 x 262 x 117 mm	2132 gm
21-465-03	580 x 280 x 150 mm	550 x 262 x 132 mm	2270 gm
21-465-04	580 x 280 x 200 mm	550 x 262 x 182 mm	2627 gm
21-465-05	580 x 280 x 260 mm	550 x 262 x 244 mm	2851 gm




Perforated lid



Grey	21-466-01
Yellow	21-466-02
Red	21-466-03
Blue	21-466-04
Green	21-466-05
Black	21-466-06

Weight: 1111 gm
 Size: 580 x 280 mm

Safety Lid



Grey	21-467-01
Yellow	21-467-02
Red	21-467-03
Blue	21-467-04
Green	21-467-05
Black	21-467-06

Weight: 851.5 gm
 Size: 580 x 280 mm

1/1 Size Container - Wire Baskets



Wire Basket

21-548-01	485 x 255 x 30 mm
21-548-02	485 x 255 x 50 mm
21-548-03	485 x 255 x 70 mm
21-548-04	485 x 255 x 100 mm



Wire Basket

21-545-01	540 x 255 x 30 mm
21-545-02	540 x 255 x 50 mm
21-545-03	540 x 255 x 70 mm
21-545-04	540 x 255 x 100 mm



Wire Base

21-611-03	545 x 255 x 30 mm
-----------	-------------------



Drape Retainer

21-612-05	545 x 255 mm
-----------	--------------



Drape

21-610-07	130 x 60 cm
21-610-08	140 x 100 cm
21-610-09	185 x 45 cm
21-610-10	185 x 60 cm



Non perforated
Bottom

21-468-01
21-468-02
21-468-03
21-468-04
21-468-05
21-468-06

Weight: 2268.5 gm



Perforated
Bottom

21-469-01
21-469-02
21-469-03
21-469-04
21-469-05
21-469-06

Weight: 2503 gm

1/1 Size Container - Safety Model



Outside: 465 x 280 x 100 mm
 Inside: 425 x 265 x 84 mm

21-470-01
21-470-02
21-470-03
21-470-04
21-470-05
21-470-06

Weight: 2549.5 gm

21-471-01
21-471-02
21-471-03
21-471-04
21-471-05
21-471-06

Weight: 2715 gm



Outside: 465 x 280 x 135 mm
 Inside: 425 x 265 x 115 mm

21-472-01
21-472-02
21-472-03
21-472-04
21-472-05
21-472-06

Weight: 2620.5 gm

21-473-01
21-473-02
21-473-03
21-473-04
21-473-05
21-473-06

Weight: 2718 gm



Outside: 465 x 280 x 150 mm
 Inside: 425 x 260 x 130 mm



Non perforated
Bottom

Grey	21-474-01
Yellow	21-474-02
Red	21-474-03
Blue	21-474-04
Green	21-474-05
Black	21-474-06

Weight: 2980 gm



Perforated
Bottom

Grey	21-475-01
Yellow	21-475-02
Red	21-475-03
Blue	21-475-04
Green	21-475-05
Black	21-475-06

Weight: 3213.5 gm



Outside: 465 x 280 x 100 mm
Inside: 425 x 265 x 84 mm

1/1 Size Container - Safety Model

Grey	21-476-01
Yellow	21-476-02
Red	21-476-03
Blue	21-476-04
Green	21-476-05
Black	21-476-06

Weight: 3260 gm

Grey	21-477-01
Yellow	21-477-02
Red	21-477-03
Blue	21-477-04
Green	21-477-05
Black	21-477-06

Weight: 2425 gm



Outside: 465 x 280 x 135 mm
Inside: 425 x 265 x 115 mm

Grey	21-478-01
Yellow	21-478-02
Red	21-478-03
Blue	21-478-04
Green	21-478-05
Black	21-478-06

Weight: 2975 gm

Grey	21-479-01
Yellow	21-479-02
Red	21-479-03
Blue	21-479-04
Green	21-479-05
Black	21-479-06

Weight: 3561 gm



Outside: 465 x 280 x 150 mm
Inside: 425 x 260 x 130 mm



Non perforated
Bottom

	Outside	Inside	Weight
21-480-01	465 x 280 x 100 mm	465 x 262 x 84 mm	1396 gm
21-480-02	465 x 280 x 135 mm	465 x 262 x 117 mm	1676 gm
21-480-03	465 x 280 x 150 mm	465 x 262 x 132 mm	1792 gm

1/1 Size Container - Bottoms and Lids



Perforated
Bottom

	Outside	Inside	Weight
21-481-01	465 x 280 x 100 mm	425 x 262 x 84 mm	2044 gm
21-481-02	465 x 280 x 135 mm	425 x 262 x 117 mm	2132 gm
21-481-03	465 x 280 x 150 mm	425 x 262 x 132 mm	2270 gm




Perforated lid



21-482-01
21-482-02
21-482-03
21-482-04
21-482-05
21-482-06

Weight: 873 gm
 Size: 465 x 280 mm

Safety Lid



21-483-01
21-483-02
21-483-03
21-483-04
21-483-05
21-483-06

Weight: 710 gm
 Size: 465 x 280 mm

3/4 Size Container - Wire Baskets



Wire Basket

21-551-01	405 x 255 x 30 mm
21-551-02	405 x 255 x 50 mm
21-551-03	405 x 255 x 70 mm
21-551-04	405 x 255 x 100 mm



Wire Base

21-611-02	405 x 255 x 30 mm
-----------	-------------------



Drape Retainer


21-612-04	405 x 255 mm
-----------	--------------



Drape

21-610-08	140 x 100 cm
-----------	--------------

1/2 Size Container - Standard Model



Non perforated
Bottom



Perforated
Bottom

	21-484-01
	21-484-02
	21-484-03
	21-484-04
	21-484-05
	21-484-06

Weight: 1744.5 gm

	21-485-01
	21-485-02
	21-485-03
	21-485-04
	21-485-05
	21-485-06

Weight: 1890.5 gm



Outside: 285 x 280 x 100 mm
Inside: 262 x 260 x 84 mm

	21-486-01
	21-486-02
	21-486-03
	21-486-04
	21-486-05
	21-486-06

Weight: 1947 gm

	21-487-01
	21-487-02
	21-487-03
	21-487-04
	21-487-05
	21-487-06

Weight: 2125 gm



Outside: 285 x 280 x 135 mm
Inside: 262 x 260 x 115 mm

	21-488-01
	21-488-02
	21-488-03
	21-488-04
	21-488-05
	21-488-06

Weight: 1981 gm

	21-489-01
	21-489-02
	21-489-03
	21-489-04
	21-489-05
	21-489-06

Weight: 2462.5 gm



Outside: 285 x 280 x 150 mm
Inside: 262 x 260 x 132 mm

	21-490-01
	21-490-02
	21-490-03
	21-490-04
	21-490-05
	21-490-06

Weight: 2101 gm

	21-491-01
	21-491-02
	21-491-03
	21-491-04
	21-491-05
	21-491-06

Weight: 2513.5 gm



Outside: 285 x 280 x 200 mm
Inside: 262 x 260 x 180 mm

	21-492-01
	21-492-02
	21-492-03
	21-492-04
	21-492-05
	21-492-06

Weight: 2483.5 gm

	21-493-01
	21-493-02
	21-493-03
	21-493-04
	21-493-05
	21-493-06

Weight: 2657.5 gm



Outside: 285 x 280 x 260 mm
Inside: 262 x 260 x 238 mm

1/2 Size Container - Safety Model



Non perforated
Bottom



Perforated
Bottom

21-494-01
21-494-02
21-494-03
21-494-04
21-494-05
21-494-06

Weight: 2162.5 gm

21-495-01
21-495-02
21-495-03
21-495-04
21-495-05
21-495-06

Weight: 2309.5 gm



Outside: 285 x 280 x 100 mm
Inside: 262 x 260 x 84 mm

21-496-01
21-496-02
21-496-03
21-496-04
21-496-05
21-496-06

Weight: 2351 gm

21-497-01
21-497-02
21-497-03
21-497-04
21-497-05
21-497-06

Weight: 2466 gm



Outside: 285 x 280 x 135 mm
Inside: 262 x 260 x 116 mm

21-498-01
21-498-02
21-498-03
21-498-04
21-498-05
21-498-06

Weight: 2394 gm

21-499-01
21-499-02
21-499-03
21-499-04
21-499-05
21-499-06

Weight: 2583 gm



Outside: 285 x 280 x 150 mm
Inside: 262 x 260 x 132 mm

21-500-01
21-500-02
21-500-03
21-500-04
21-500-05
21-500-06

Weight: 2560 gm

21-501-01
21-501-02
21-501-03
21-501-04
21-501-05
21-501-06

Weight: 2785.5 gm



Outside: 285 x 280 x 200 mm
Inside: 262 x 260 x 182 mm

21-502-01
21-502-02
21-502-03
21-502-04
21-502-05
21-502-06

Weight: 2930 gm

21-503-01
21-503-02
21-503-03
21-503-04
21-503-05
21-503-06

Weight: 3069.5 gm



Outside: 285 x 280 x 260 mm
Inside: 262 x 260 x 238 mm



Non perforated
Bottom

	Outside	Inside	Weight
21-504-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1096 gm
21-504-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1284 gm
21-504-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1327 gm
21-504-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1447 gm
21-504-05	285 x 280 x 260 mm	262 x 262 x 240 mm	1817 gm

1/2 Size Container - Bottoms and Lids



Perforated
Bottom

	Outside	Inside	Weight
21-505-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1242 gm
21-505-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1398 gm
21-505-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1471 gm
21-505-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1672 gm
21-505-05	285 x 280 x 260 mm	262 x 262 x 240 mm	2000 gm



Perforated lid



21-506-01
21-506-02
21-506-03
21-506-04
21-506-05
21-506-06

Weight: 254 gm
 Size: 285 x 280 mm

Safety Lid



21-507-01
21-507-02
21-507-03
21-507-04
21-507-05
21-507-06

Weight: 458 gm
 Size: 285 x 280 mm

1/2 Size Container - Wire Baskets



Wire Basket

21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm
21-554-03	255 x 245 x 70 mm
21-554-04	255 x 245 x 100 mm



Wire Base

21-611-01	255 x 245 x 30 mm
-----------	-------------------



Drape Retainer

21-612-03	255 x 245 mm
-----------	--------------



Drape

21-610-04	90 x 60 cm
21-610-05	90 x 80 cm
21-610-06	130 x 45 cm
21-610-07	130 x 60 cm

Flat Container - Containers and Bottoms



Non perforated
Bottom



Perforated
Bottom

	21-508-01
	21-508-02
	21-508-03
	21-508-04
	21-508-05
	21-508-06

Weight: 1285 gm



Outside: 285 x 280 x 55 mm
 Inside: 270 x 270 x 33 mm

	21-509-01
	21-509-02
	21-509-03
	21-509-04
	21-509-05
	21-509-06

Weight: 1377 gm

	21-510-01
	21-510-02
	21-510-03
	21-510-04
	21-510-05
	21-510-06

Weight: 1529 gm



Outside: 600 x 400 x 180 mm
 Inside: 550 x 350 x 160 mm



Non perforated
Bottom

	Outside	Inside	Weight
21-511-01	285 x 280 x 55 mm	270 x 270 x 33 mm	611 gm
21-511-02	285 x 280 x 85 mm	270 x 270 x 64 mm	703 gm



Perforated
Bottom

	Outside	Inside	Weight
21-512-01	285 x 280 x 85 mm	270 x 270 x 64 mm	857 gm



Flat Container - Lids and accessories

Perforated lid



21-513-01
21-513-02
21-513-03
21-513-04
21-513-05
21-513-06

Weight: 674 gm



Wire Basket

21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm




Drape Retainer

21-612-03	255 x 245 mm
-----------	--------------



Wire Base

21-611-01	255 x 245 x 30 mm
-----------	-------------------



Silicone Mat

21-628-02	245 x 245 mm
-----------	--------------



Drape

21-610-04	90 x 60 cm
21-610-05	90 x 80 cm

Mini Container



Non perforated
Bottom

Grey	21-514-01
Yellow	21-514-02
Red	21-514-03
Blue	21-514-04
Green	21-514-05
Black	21-514-06

Weight: 510.5 gm



Perforated
Bottom

Grey	21-515-01
Yellow	21-515-02
Red	21-515-03
Blue	21-515-04
Green	21-515-05
Black	21-515-06

Weight: 589 gm



Outside: 300 x 140 x 40 mm
 Inside: 280 x 130 x 27 mm

Grey	21-516-01
Yellow	21-516-02
Red	21-516-03
Blue	21-516-04
Green	21-516-05
Black	21-516-06

Weight: 567 gm

Grey	21-517-01
Yellow	21-517-02
Red	21-517-03
Blue	21-517-04
Green	21-517-05
Black	21-517-06

Weight: 630 gm



Outside: 300 x 140 x 70 mm
 Inside: 280 x 130 x 55 mm

Grey	21-518-01
Yellow	21-518-02
Red	21-518-03
Blue	21-518-04
Green	21-518-05
Black	21-518-06

Weight: 625.5 gm

Grey	21-519-01
Yellow	21-519-02
Red	21-519-03
Blue	21-519-04
Green	21-519-05
Black	21-519-06

Weight: 625 gm



Outside: 300 x 140 x 100 mm
 Inside: 280 x 130 x 84 mm



Non perforated
Bottom


	Outside	Inside	Weight
21-520-01	300 x 140 x 40 mm	280 x 130 x 27 mm	190 gm
21-520-02	300 x 140 x 70 mm	280 x 130 x 55 mm	241 gm
21-520-03	300 x 140 x 100 mm	280 x 130 x 84 mm	305 gm



Perforated
Bottom

	Outside	Inside	Weight
21-521-01	300 x 140 x 40 mm	280 x 130 x 27 mm	258 gm
21-521-02	300 x 140 x 70 mm	280 x 130 x 55 mm	304 gm
21-521-03	300 x 140 x 100 mm	280 x 130 x 84 mm	326 gm

Perforated lid



	21-522-01
	21-522-02
	21-522-03
	21-522-04
	21-522-05
	21-522-06

Weight: 325.5 gm

Mini Container - Bottoms and Lids



Mini Container - Accessories



Drape Retainer

21-612-01 280 x 125 mm



Silicone Mat

21-628-01 270 x 125 mm



Drape

21-610-01 60 x 40 cm

Dental Containers



Perforated
Bottom

21-523-01

Weight: 783 gm



Outside: 310 x 190 x 40 mm
Inside: 280 x 182 x 35 mm



Perforated
Bottom

21-524-01

Weight: 862 gm



Outside: 310 x 190 x 65 mm
Inside: 280 x 182 x 55 mm



Perforated
Bottom

21-525-01

Weight: 1022 gm



Outside: 310 x 190 x 130 mm
Inside: 280 x 182 x 124 mm

Dental Container - Bottom, Lid and Accessories




Perforated
Bottom

	Outside	Inside	Weight
21-526-01	310 x 190 x 40 mm	280 x 182 x 35 mm	341 gm
21-526-02	310 x 190 x 65 mm	280 x 182 x 55 mm	415 gm
21-526-03	310 x 190 x 130 mm	280 x 182 x 124 mm	570 gm



Perforated lid



21-527-01

Weight: 417 gm
 Size: 310 x 190 mm



Drape Retainer

21-612-02 310 x 190 mm



Drape

21-610-01	60 x 40 cm
21-610-02	65 x 65 cm
21-610-03	85 x 85 cm

Implant Containers



Lid perforated
Bottom perforated (no filter)

	21-528-01
	21-528-02
	21-528-03
	21-528-04
	21-528-05
	21-528-06

Weight: 654 gm



Outside: 500 x 155 x 75 mm
Inside: 495 x 155 x 75 mm



Lid perforated
Bottom perforated (no filter)

	21-529-01
	21-529-02
	21-529-03
	21-529-04
	21-529-05
	21-529-06

Weight: 368 gm



Outside: 500 x 155 x 75 mm
Inside: 495 x 155 x 75 mm

Optik Containers



Lid perforated
 Bottom perforated (no filter)

21-530-01

Weight: 302 gm



Outside: 450 x 70 x 70 mm
 Inside: 450 x 70 x 30 mm



Lid perforated
 Bottom perforated (no filter)


21-531-01

Weight: 152 gm



Outside: 250 x 60 x 50 mm
 Inside: 250 x 55 x 50 mm

Removal Containers




Non perforated
Lid and Bottom

	21-532-01
	21-532-02
	21-532-03
	21-532-04
	21-532-05
	21-532-06

Weight: 1677 gm



Outside: 285 x 280 x 100 mm
Inside: 262 x 260 x 82 mm




Non perforated
Lid and Bottom

	21-533-01
	21-533-02
	21-533-03
	21-533-04
	21-533-05
	21-533-06

Weight: 1929 gm



Outside: 285 x 280 x 150 mm
Inside: 262 x 260 x 132 mm



Non perforated
Lid and Bottom

	21-534-01
	21-534-02
	21-534-03
	21-534-04
	21-534-05
	21-534-06

Weight: 2883 gm



Outside: 580 x 280 x 150 mm
Inside: 548 x 258 x 130 mm

Removal Container - Bottoms and Lids



Non perforated
Bottom

	Outside	Weight
21-535-01	580 x 280 x 150 mm	2108 gm



Non perforated
Bottom

	Outside	Weight
21-536-01	285 x 280 x 150 mm	1096 gm
21-636-02	285 x 280 x 150 mm	1327 gm




Non-Perforated lid



	21-537-01
	21-537-02
	21-537-03
	21-537-04
	21-537-05
	21-537-06

Weight: 458 gm
 Size: 285 x 280 mm

Non-Perforated lid



	21-538-01
	21-538-02
	21-538-03
	21-538-04
	21-538-05
	21-538-06

Weight: 875 gm
 Size: 580 x 280 mm

Brush Container - Stainless Steel and Aluminium



21-539-10 For 10 brushes

21-539-20 For 20 brushes

Stainless steel

	21-540-01
	21-540-02
	21-540-03
	21-540-04
	21-540-05

Aluminium
For 20 brushes
Size: 240 x 275 x 70 mm



Single Tray Containers

21-541-01	200 x 140 x 50 mm
21-541-02	280 x 195 x 50 mm
21-541-03	265 x 255 x 50 mm
21-541-04	380 x 260 x 50 mm

For micro instruments



Double Tray Containers

21-542-01	280 x 195 x 50 mm
21-542-02	265 x 255 x 50 mm
21-542-03	380 x 260 x 50 mm

With 2 take-out perforated trays

Baskets



21-543-01	540 x 255 x 30 mm
21-543-02	540 x 255 x 50 mm
21-543-03	540 x 255 x 70 mm
21-543-04	540 x 255 x 100 mm

Wire basket with lid

21-544-00

Lid only



21-545-01	540 x 255 x 30 mm
21-545-02	540 x 255 x 50 mm
21-545-03	540 x 255 x 70 mm
21-545-04	540 x 255 x 100 mm

Wire basket



21-546-01	485 x 255 x 30 mm
21-546-02	485 x 255 x 50 mm
21-546-03	485 x 255 x 70 mm
21-546-04	485 x 255 x 100 mm

Wire basket with lid

21-547-00

Lid only



21-548-01	485 x 255 x 30 mm
21-548-02	485 x 255 x 50 mm
21-548-03	485 x 255 x 70 mm
21-548-04	485 x 255 x 100 mm

Wire basket

Baskets



21-549-01	405 x 255 x 30 mm
21-549-02	405 x 255 x 50 mm
21-549-03	405 x 255 x 70 mm
21-549-04	405 x 255 x 100 mm

Wire basket with lid

21-550-00

Lid only



21-551-01	405 x 255 x 30 mm
21-551-02	405 x 255 x 50 mm
21-551-03	405 x 255 x 70 mm
21-551-04	405 x 255 x 100 mm

Wire basket



21-552-01	255 x 245 x 30 mm
21-552-02	255 x 245 x 50 mm
21-552-03	255 x 245 x 70 mm
21-552-04	255 x 245 x 100 mm

Wire basket with lid

21-553-00


Lid only



21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm
21-554-03	255 x 245 x 70 mm
21-554-04	255 x 245 x 100 mm

Wire basket

Baskets



21-555-01 540 x 255 x 50 mm

21-555-02 540 x 255 x 70 mm

21-555-03 540 x 255 x 100 mm

Wire basket with lid

21-556-00

Lid only

21-557-01 485 x 255 x 50 mm

21-557-02 485 x 255 x 70 mm

21-557-03 485 x 255 x 100 mm

Wire basket with lid

21-558-00

Lid only



21-561-01 405 x 255 x 50 mm

21-561-02 405 x 255 x 70 mm

21-561-03 405 x 255 x 100 mm

Wire basket with lid

21-562-00

Lid only



21-564-01 255 x 245 x 50 mm


21-564-02 255 x 245 x 70 mm

21-564-03 255 x 245 x 100 mm

Wire basket with lid

21-565-00

Lid only



21-559-01 540 x 255 x 50 mm

21-559-02 540 x 255 x 70 mm

21-559-03 540 x 255 x 100 mm

Wire basket

21-560-01 485 x 255 x 50 mm

21-560-02 485 x 255 x 70 mm

21-560-03 485 x 255 x 100 mm

Wire basket



21-563-01 405 x 255 x 50 mm

21-563-02 405 x 255 x 70 mm

21-563-03 405 x 255 x 100 mm

Wire basket



21-566-01 255 x 245 x 50 mm

21-566-02 255 x 245 x 70 mm

21-566-03 255 x 245 x 100 mm

Wire basket

Baskets



21-567-01 290 x 80 x 52 mm

21-567-02 460 x 80 x 52 mm

21-567-03 670 x 80 x 52 mm

Endoscope basket



21-568-00 640 x 150 x 77 mm

Endoscope basket



21-569-00 274 x 172 x 43 mm

Wire basket



21-570-00 267 x 125 x 50 mm

Wire basket



21-571-00 274 x 172 x 43 mm


Wire basket



21-572-00 267 x 125 x 50 mm

Wire basket

Wire Baskets



21-573-01	30 x 30 x 30 mm
21-573-02	40 x 40 x 20 mm
21-573-03	40 x 40 x 40 mm
21-573-04	60 x 60 x 20 mm
21-573-05	60 x 60 x 40 mm
21-573-06	60 x 60 x 60 mm
21-573-07	80 x 80 x 20 mm
21-573-08	80 x 40 x 40 mm
21-573-09	80 x 80 x 40 mm
21-573-10	80 x 80 x 60 mm
21-573-11	80 x 80 x 80 mm
21-573-12	100 x 100 x 10 mm
21-573-13	120 x 120 x 30 mm
21-573-14	120 x 120 x 60 mm
21-573-15	120 x 120 x 80 mm
21-573-16	240 x 240 x 30 mm
21-573-17	240 x 240 x 60 mm
21-573-18	240 x 240 x 80 mm
21-573-19	240 x 240 x 100 mm

Small fine mesh basket with lid
Mesh distance, 2 x 2 mm



21-574-01	600 x 400 x 150 mm
21-574-02	600 x 400 x 300 mm

Laundry basket



21-575-01	600 x 400 x 150 mm
21-575-02	600 x 400 x 300 mm

Laundry basket

Baskets



Size: 1/1

Box and lid	Box only	Lid only	
21-576-01	21-577-01	21-578-01	540 x 253 x 30 mm
21-576-02	21-577-02	21-578-02	540 x 253 x 50 mm
21-576-03	21-577-03	21-578-03	540 x 253 x 70 mm
21-576-04	21-577-04	21-578-04	540 x 253 x 100 mm

Size: 1/1

Box and lid	Box only	Lid only	
21-579-01	21-580-01	21-581-01	485 x 253 x 30 mm
21-579-02	21-580-02	21-581-02	485 x 253 x 50 mm
21-579-03	21-580-03	21-581-03	485 x 253 x 70 mm
21-579-04	21-580-04	21-581-04	485 x 253 x 100 mm

Size: 3/4

Box and lid	Box only	Lid only	
21-582-01	21-583-01	21-584-01	405 x 253 x 30 mm
21-582-02	21-583-02	21-584-02	405 x 253 x 50 mm
21-582-03	21-583-03	21-584-03	405 x 253 x 70 mm
21-582-04	21-583-04	21-584-04	405 x 253 x 100 mm

Size: 1/2

Box and lid	Box only	Lid only	
21-585-01	21-586-01	21-587-01	244 x 253 x 30 mm
21-585-02	21-586-02	21-587-02	244 x 253 x 50 mm
21-585-03	21-586-03	21-587-03	244 x 253 x 70 mm
21-585-04	21-586-04	21-587-04	244 x 253 x 100 mm

Size: Dental

Box and lid	Box only	Lid only	
21-588-01	21-589-01	21-590-01	275 x 175 x 40 mm
21-588-02	21-589-02	21-590-02	275 x 175 x 95 mm

Size: Small

Box and lid	Box only	Lid only	
21-591-01	21-592-01	21-593-01	265 x 125 x 50 mm




21-594-14
 Stand
 14 mm height




21-595-20
 Stand
 20 mm height

Filters



21-596-00 210 x 95 mm

Long term textile filter for mini and dental container for 50 sterilization cycles.



21-597-00 235 x 118 mm

Long term PTFE filter for mini and dental container for 1200 sterilization cycles.




21-598-00 235 x 118 mm

Long term PTFE filter for 1/1, 1/2, 3/4, mini and dental container.



21-599-00 Ø 190 mm

Long term textil filter for 1/1, 1/2, 3/4 container for 50 sterilization cycles.



21-600-00 235 x 118 mm

Long term PTFE filter for 1/1, 1/2, 3/4 container for 1200 sterilization cycles.



21-601-00 235 x 118 mm

Disposable filter for 1/1, 1/2, 3/4 container.

Indicators, Labels and Drapes



21-602-00

Paper label with indicator
 For 1/1, 1/2, 3/4 containers
 Pack of 100



21-602-00

Paper label with indicator
 For mini and dental containers
 Pack of 100



21-602-00

Heavy duty clip
 Serrated jaws



21-605-01 For paper filter

21-605-02 For textile filter

Filter tray



21-606-00

Filter dispenser



21-607-00

Plastic security seal
 Pack of 100




Drape

21-610-01	60 x 40 cm
21-610-02	65 x 65 cm
21-610-03	85 x 85 cm
21-610-04	90 x 60 cm
21-610-05	90 x 80 cm
21-610-06	130 x 45 cm
21-610-07	130 x 60 cm
21-610-08	140 x 100 cm
21-610-09	185 x 45 cm
21-610-10	185 x 60 cm



21-607-00

Dispenser for
 identification labels



Coloured identification labels.
 Printed to your specification,
 16 letters max.

Wire Base, Drape Retainer, Metal Labels and Fixation Clamps



Wire Base

21-611-01	255 x 245 x 30 mm
21-611-02	405 x 255 x 30 mm
21-611-03	545 x 255 x 30 mm



Drape Retainer

21-612-01	280 x 125 mm
21-612-02	310 x 190 mm
21-612-03	255 x 245 mm
21-612-04	405 x 255 mm
21-612-05	545 x 255 mm



Metal Label for Wire Basket

21-613-00	545 x 255 x 30 mm
-----------	-------------------

Unlabeled



Metal Label for Wire Basket

21-614-00	545 x 255 mm
-----------	--------------

Labeled



21-615-01

24 x 28 x 25 mm
Ø 7 - 12 mm



21-615-02

24 x 28 x 29 mm
Ø 12 - 16 mm



21-615-03

24 x 28 x 31 mm
Ø 16 - 26 mm



21-615-04

24 x 28 x 40 mm
Ø 26 - 36 mm



21-615-05

24 x 28 x 44 mm
Ø 36 - 45 mm



21-616-01

50 x 20 x 12 mm



21-616-02

65 x 20 x 20 mm



Tray Tags



	21-617-01
	21-617-02
	21-617-03
	21-617-04
	21-617-05

	21-617-01
	21-617-02
	21-617-03
	21-617-04
	21-617-05

Nylon
 Size: 105x40mm
 Pack of 100



21-618-00
 Stainless steel
 Attachment ring
 Ø 30 mm
 Pack of 100

Components




21-619-00

Holding pin for
partition sheets



21-620-00

Holding pin for
wire baskets



21-621-01 50 x 20 mm

21-621-02 130 x 20 mm

21-621-03 225 x 20 mm

21-621-04 460 x 20 mm

Sheets



21-622-01 230 mm

21-622-02 470 mm

Sliding rails with securing screws



21-623-01 230 mm

21-623-02 470 mm

Limiting bar with securing screws



21-624-00

Divider, to be used
with 21-622-01 and 21-622-02



21-625-00


Spiral holding



21-626-00

Double spiral holding

Silicone Mats




Silicone Mat Fig. 1

21-627-01	270 x 125 mm
21-627-02	245 x 245 mm
21-627-03	220 x 230 mm
21-627-04	380 x 230 mm
21-627-05	440 x 230 mm
21-627-06	520 x 230 mm



Silicone Mat Fig. 2

21-628-01	270 x 125 mm
21-628-02	245 x 245 mm
21-628-03	220 x 230 mm
21-628-04	380 x 230 mm
21-628-05	440 x 230 mm
21-628-06	520 x 230 mm



Silicone Mat Fig. 3

21-629-00	545 x 255 mm
-----------	--------------





Silicone Mat Fig. 4

21-630-01	220 x 125 mm
21-630-02	380 x 230 mm
21-630-03	440 x 230 mm
21-630-04	520 x 230 mm


Silicone Fixation System





	Length	Height
21-631-01	97 mm	41 mm
21-631-02	125 mm	41 mm
21-631-03	152 mm	41 mm
21-631-04	250 mm	41 mm




	Length	Height
21-632-01	90 mm	41 mm
21-632-02	120 mm	41 mm
21-632-03	150 mm	41 mm
21-632-04	240 mm	41 mm



	Length	Height
21-633-01	87 mm	41 mm
21-633-02	119 mm	41 mm
21-633-03	150 mm	41 mm
21-633-04	238 mm	41 mm




	Length	Height
21-634-01	77 mm	41 mm
21-634-02	112 mm	41 mm
21-634-03	158 mm	41 mm
21-634-04	238 mm	41 mm




	Length	Height
21-635-00	178 mm	20 mm

Instruments holder from Ø 3 to 8 mm




	Length	Height
21-636-00	73 mm	20 mm

Instruments holder from Ø 5




	Length	Height
21-637-00	73 mm	20 mm

Instruments holder from Ø 5



	Length	Height
21-638-00	178 mm	30 mm



Instruments holder from Ø 3 to 8 mm




Silicone Fixation System




Length Height
 21-639-00 121 mm 18 mm




Length Height
 21-640-00 123 mm 18 mm




Length Height
 21-641-00 178 mm 20 mm
 Instruments holder from Ø 3 mm





Length Height
 21-642-00 123 mm 18 mm
 Instruments holder from Ø 5 mm





Length Height
 21-643-00 123 mm 18 mm





Length Height
 21-644-00 23.5 mm 21 mm
 Ø 8 mm




Length Height
 21-645-00 23.5 mm 21 mm
 Ø 5 mm




Length Height
 21-646-00 23.5 mm 21 mm
 Ø 13 mm



Length Height
 21-647-00 23.5 mm 21 mm
 Ø 10.5 mm




Silicone Fixation System




	Length	Height
21-648-00	73 mm	30 mm

Endoscope holder, double
one side Ø 15 mm open
one side from Ø 3 to Ø 15 mm



	Length	Height
21-649-00	72 mm	30 mm

Endoscope holder, double
one side Ø 15 mm open
one side from Ø 3 to Ø 15 mm




	Length	Height
21-650-00	72 mm	30 mm

Instruments holder
from Ø 3 to Ø 15 mm




	Length	Height
21-651-00	72 mm	30 mm

Instruments holder
from Ø 3 to Ø 8 mm



	Length	Height
21-652-00	38 mm	20 mm

Endoscope holder, double
one side 12 mm open
one side from Ø 3 to Ø 8 mm



	Length	Height
21-638-00	72 mm	30 mm


Endoscope holder, double
one side 12 mm open
one side from Ø 3 to Ø 8 mm

Silicone Fixation System




	Length	Height
21-654-00	54 mm	20 mm

Endoscope holder
 Ø 4 mm closed




	Length	Height
21-655-00	38 mm	29 mm

Endoscope holder
 Ø 6 mm closed




	Length	Height
21-656-00	38 mm	29 mm

Endoscope holder
 Ø 8 mm closed




	Length	Height
21-657-00	38 mm	29 mm

Endoscope holder
 Ø 12 mm closed



	Length	Height
21-658-00	38 mm	29 mm

Endoscope holder
 Ø 12 mm open



	Length	Height
21-659-00	38 mm	29 mm

Endoscope holder
 Ø 15 mm open